

STORZ

KARL STORZ — ENDOSKOPE

en **Instructions for use**
**Uretero-renoscope, ureteroscope, cystoscope-
urethroscope**



10-2020

Copyright ©

All product illustrations, product descriptions, and texts are the intellectual property of KARL STORZ SE & Co. KG.

Their use and reproduction by third parties require the express approval of KARL STORZ SE & Co. KG.

All rights reserved.

Table of contents

1	General information.....	4
1.1	Read the instructions for use.....	4
1.2	Read the instructions for use of combinable products.....	4
1.3	Scope.....	4
1.4	Description of warning messages.....	4
2	Normal use.....	6
2.1	Intended use	6
2.2	Indications.....	6
2.3	Contraindications.....	8
2.4	Target user populations	8
2.5	Patient groups.....	9
3	Safety.....	10
3.1	Serious incidents	10
3.2	Unsterile product	10
3.3	Correct handling	10
3.4	Damaged products	10
3.5	Working in the field of vision.....	10
3.6	Hot components	11
3.7	High light intensity	11
3.8	Appropriate combination of endoscope and fiber optic light cable	11
3.9	Patient leakage current.....	11
3.10	Creutzfeldt-Jakob disease	11
4	Product description	12
4.1	Product overview	12
4.2	Possible combinations.....	14
4.3	Ambient conditions.....	16
4.4	Symbols on the packaging	17
5	Preparation.....	18
5.1	Unpacking the product.....	18
5.2	Assembling the product.....	18
6	Disassembly.....	19
6.1	Disassembling the product.....	19
7	Maintenance, servicing, repairs, and disposal.....	20
7.1	Repairing the product	20
7.2	Disposing of the product	20
8	Accessories and spare parts	21
8.1	Accessories.....	21
8.2	Spare parts	21
9	Subsidiaries.....	22

1 General information

1.1 Read the instructions for use

If the instructions for use are not followed, patients, users, and third parties may be injured or the product may be damaged.

- ▶ Read the instructions for use carefully and follow all the safety notes and warnings.
- ▶ Keep the instructions for use clearly visible next to the product.

1.2 Read the instructions for use of combinable products

If the instructions for use of combinable products are not followed, patients, users, and third parties may be injured or the product may be damaged.

- ▶ Read the instructions for use of the combinable products carefully and follow all the safety notes and warnings.

1.3 Scope

These instructions for use are valid for the following products:

Product	Item number
Ureteroscope	27000K
	27001K
	27002K
	27011K
Uretero-roscope	27000L
	27001L
	27002L
	27003L
	27011L
	27013L
Uretero-roscope, pediatric	27002KP
Cystoscope-urethroscope, pediatric	27030KA/KB
Instrument inlet	27001G, 27001GF 27001GH
Insertion aid for guide wires	27001E

1.4 Description of warning messages

To prevent any injury to persons or damage to property, the warnings and safety notes in the instructions for use must be observed. The warning messages describe the following levels of danger.

▲ WARNING

WARNING

Designates a possible imminent risk. If this is not avoided, it could lead to death or serious injuries.

▲ CAUTION

CAUTION

Designates a possible imminent risk. If this is not avoided, it could lead to minor injuries.

ATTENTION

ATTENTION

Designates a possibly harmful situation. If this is not avoided, the products could be damaged.

2 Normal use

2.1 Intended use

Ureterscopes, uretero-renoscopes (27000K, 27001K, 27002K, 27011K, 27000L, 27001L, 27002L, 27003L, 27011L, 27013L)

Rigid/semirigid ureterscopes/uretero-renoscopes (telescope) are used for endoscopic imaging of the ureter, the pyeloureteral junction, and the renal pelvis for the purposes of diagnosis and treatment.

Ureterscopes/uretero-renoscopes with an integrated working channel allow for the use of auxiliary instruments such as forceps, scissors, stone baskets, etc., laser fibers and other lithotripters for the treatment of urinary stones, stenoses/strictures, and benign/malignant tumors. Ureterscopes/uretero-renoscopes are designed for short-term use in invasive procedures through a body orifice.

Pediatric uretero-renoscopes (27002KP)

The uretero-renaloscope/ureterscope (telescope) or flex. uretero-renaloscope is used for endoscopic imaging of the ureter and renal collecting system for the purposes of endourological diagnosis and treatment. The instruments are designed for short-term use in invasive procedures through a body orifice.

Pediatric cystoscope-urethroscope (27030KA/KB)

The telescope is used for endoscopic imaging and, if necessary, for creating a working or irrigation channel during endourological diagnosis and treatment. Compact cystoscopes are used for endoscopic imaging of the urethra and bladder for the purposes of diagnosis and treatment. The integrated working channel makes it possible to use flexible or semirigid auxiliary instruments such as forceps, stone baskets, etc., coagulation electrodes and laser fibers for lithotripsy or for the ablation of benign/malignant tissue. Telescopes/compact cystoscopes are designed for transient use in invasive procedures through a body orifice.

2.2 Indications

Ureterscopes, uretero-renoscopes (27000K, 27001K, 27002K, 27011K, 27000L, 27001L, 27002L, 27003L, 27011L, 27013L)

The use of products for ureteroscopy/ureterorenoscopy is indicated if uretero(reno)scopy is indicated in the opinion of the responsible physician.

This includes:

- Further diagnostic measures in the event of conspicuous imaging and filling defects in the upper urinary tract shown in the imaging
- Pathological bladder-wash cytology without reference to bladder tumors or from the upper urinary tract
- Unclear hematuria (macro/microhematuria)
- Diagnosis and treatment of urinary stones in the ureter or renal collecting system
- Diagnosis and treatment of ureteral strictures or tumors – iatrogenically introduced foreign material
- Narrow pyeloureteral junction
- Possibly the percutaneous antegrade treatment of kidney stones, tumors, and ureter outflow constrictions
- UUT (Upper Urinary Tract) tumors or UTUC (Upper Tract Urothelial Cancer)

- Nephrolithiasis of lower calyceal group
- Positive bladder-wash cytology from the upper urinary tract
- Tumors in the upper urinary tract

Indication for retrograde treatment of a transitional cell carcinoma:

- Impaired renal function/single kidney
- Bilateral tumor development
- Multiple comorbidities/limited life expectancy
- Tumor stages in which only palliation is possible

Pediatric uretero-renoscopes (27002KP)

The use of instruments for transurethral/percutaneous pediatric urology is indicated if an endourological procedure is indicated in the opinion of the responsible physician.

This includes, among other things:

Cystoscopy:

- Anomalies in the ureteral orifices
- Ureteroceles
- Vesicorenal reflux
- Enuresis
- Neurogenic micturition disorder/reflexive bladder

Ureterorenoscopy/PCNL:

- Nephrolithiasis/urolithiasis
- Ureteral strictures/stenoses
- Double ureter
- Upper ureteral calculi

Pediatric cystoscope-urethroscope (27030KA/KB)

The use of instruments for transurethral/percutaneous pediatric urology is indicated if an endourological procedure is indicated in the opinion of the responsible physician.

This includes, among other things:

Cystoscopy

- Anomalies in the ureteral orifices
- Ureteroceles
- Vesicorenal reflux
- Enuresis
- Neurogenic micturition disorder/reflexive bladder

Ureterorenoscopy/PCNL:

- Nephrolithiasis/urolithiasis
- Ureteral strictures/stenoses
- Double ureter
- Upper ureteral calculi

2.3 Contraindications

Ureteroscopes, uretero-renoscopes (27000K, 27001K, 27002K, 27011K, 27000L, 27001L, 27002L, 27003L, 27011L, 27013L)

The use of products for uretero(reno)scopy is contraindicated if, in the opinion of the responsible physician, the surgical method is contraindicated or the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Uretero(reno)scopes and their accessories must not be used for procedures in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Urinary tract infection
- Urosepsis
- Untreated/unresolved coagulation disorder
- Lithotomy position contraindicated
- Pregnancy

Pediatric uretero-renoscopes (27002KP)

The pediatric urological, transurethral, and percutaneous use of the instruments is contraindicated if, in the opinion of the responsible physician, the surgical method is contraindicated or the child is not able to undergo surgery or anesthesia due to his or her general condition. Instruments for pediatric urology must not be used for procedures in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Acute urinary tract infection
- Unresolved coagulation disorders
- Neurogenic detrusor instability

Where children are concerned, the responsible physician must question the necessity and potential therapeutic consequences of any diagnostic or therapeutic endoscopic procedures.

Pediatric cystoscope-urethroscope (27030KA/KB)

The pediatric urological, transurethral, and percutaneous use of the instruments is contraindicated if, in the opinion of the responsible physician, the surgical method is contraindicated or the child is not able to undergo surgery or anesthesia due to his or her general condition. Instruments for pediatric urology must not be used for procedures in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Acute urinary tract infection
- Unresolved coagulation disorders
- Neurogenic detrusor instability

Where children are concerned, the responsible physician must question the necessity and potential therapeutic consequences of any diagnostic or therapeutic endoscopic procedures.

2.4 Target user populations

The medical device may only be used by doctors and medical assistants with a relevant specialist qualification.

2.5 Patient groups

There are no restrictions in terms of patient groups for this product.

3 Safety

3.1 Serious incidents

According to the Medical Device Regulation (MDR), a “serious incident” includes incidents that directly or indirectly had, could have had, or could have any of the following consequences (MDR, Art. 2, No. 65 [1]):

- Death of a patient, user, or another person
- Temporary or permanent serious deterioration in the medical condition of a patient, user, or another person
- A serious threat to public health
- ▶ The manufacturer and appropriate authority must be notified of all serious incidents.

3.2 Unsterile product

The product is not sterile when delivered. The use of non-sterile products poses a risk of infection for patients, users, and third parties.

- ▶ Reprocess the product in line with the reprocessing instructions before initial use and every subsequent use.

3.3 Correct handling

If the product is not handled correctly, patients, users, and third parties may be injured.

- ▶ Only persons with the necessary medical qualification and who are acquainted with the application of the product may work with it.
- ▶ Check that the product is suitable for the procedure prior to use.
- ▶ Check the product for the following points before and after every use:
 - Completeness
 - Good working order
 - Rough surfaces left inadvertently
 - Sharp corners
 - Burred edges
 - Correct assembly of the components
 - Functionality
- ▶ Do not leave broken-off components inside the patient.
- ▶ Do not overload the product with mechanical stress.
- ▶ Do not bend bent products back to their original position.

3.4 Damaged products

Damaged products can result in injury to patients, users, or third parties.

- ▶ Before each use, check all components of the product for damage.
- ▶ Do not use damaged products.

3.5 Working in the field of vision

Using the product outside the field of vision can cause injury to tissue or can damage the product.

- ▶ Only use the product in the field of vision.

3.6 Hot components

The high level of light intensity produced by the light source may cause the distal end, the light connections, and adjacent components to heat up. This can cause burns to patients, users, and third parties.

- ▶ Set the light source output to a level that is just high enough to ensure optimal illumination of the operating area.
- ▶ Avoid contact with the distal end and light connections.

3.7 High light intensity

The high level of light intensity produced by the light source may lead to permanent eye damage or blindness, and may cause tissue and items facing the light output to heat up.

- ▶ Do not look into the light output.
- ▶ Set the light source output to a level that is just high enough to ensure optimal illumination of the operating area.
- ▶ Make sure the light output is sufficiently far away from tissue and operating accessories.

3.8 Appropriate combination of endoscope and fiber optic light cable

If the fiber diameter of the fiber optic light cable is too large, the endoscope light connection will heat up. This may cause damage to the endoscope.

If the fiber diameter of the fiber optic light cable is too small, not enough light will enter the endoscope. This may mean that the operating area is not sufficiently illuminated.

- ▶ Use the endoscope with a fiber optic light cable that has an appropriate fiber diameter.

3.9 Patient leakage current

Patient leakage currents from products may add up if powered products and powered endotherapy devices are used simultaneously. Excessively high leakage current levels may result in the patient becoming injured.

- ▶ Only use products of the same type together, particularly in the case of CF products.

3.10 Creutzfeldt-Jakob disease

Products that come into contact with the central nervous system can become contaminated by organic residue containing prions. Prions lead to infection with Creutzfeldt-Jakob disease.

If Creutzfeldt-Jakob disease has been diagnosed or is suspected:


- ▶ Dispose of the product properly and do not continue to use it.

4 Product description


4.1 Product overview



- | | | | |
|---|--------------------|---|---------------------------------|
| 1 | Light inlet socket | 4 | Mount for instrument inlet |
| 2 | Cock plug | 5 | Irrigation port |
| 3 | Spring cap | 6 | Symbol for light cable diameter |



- | | |
|---|--|
| 1 | Mount for oblique instrument inlet
(27011K/L, 27013L) |
|---|--|



- | | | | |
|---|--|---|--|
| 1 | External dimension of the instrument | 3 | Maximum diameter of instruments |
| 2 | Irrigation port inflow (additional feature
of 27003L) | 4 | Data matrix code for reading off the
serial and article numbers |



1 Instrument inlet (27001G)

2 Instrument inlet, 2 channels (27001GF)

3 Instrument inlet, 2 channels (27001GH)



1 Cleaning adaptor (27001RA)



1 Flow control stopcock (27504)

4.2 Possible combinations

Possible combinations: endoscope – instrument inlet

Uretero-renoscopes	Instrument inlet	Insertion aid
27000K	27001G	27001E (only in combination with 27001G)
27000L	27001GF	
27001K	27001GH	
27001L		
27002K		
27002L		
27003L		
27011K		
27011L		
27013L		
27030KA		
27030KB		
27002KP	27001G	

Possible combinations: endoscope – forceps, scissors, instruments

Endoscope	Forceps, scissors, instrument
27000K	27424F
27000L	27424P
27001K	27424R
27001L	27424Z
27002K	27424U
27002L	Laser fiber Lithotripsy probe
27003L	27023FM
27002KP	Laser fiber
27002K	27425F
27002L	27425P 27425R 27425Z 27425U Laser fiber Lithotripsy probe
27030KA	27071TJ 27071ZJ 27030EL 27770AA 27770A

Endoscope	Forceps, scissors, instrument
	27030M 27030N
27030KB	27095F 27095P 27095Z 27030M 27030N 27030EL 27770AA 27770A 27770B 27772AA 27772A

Possible combinations: endoscope – stone basket




Endoscope	Stone basket
27000K	27023LD
27000L	27023LF
27001K	27023TD
27001L	27023TF
27002K	27023KF
27002L	
27003L	
27011K	
27011L	
27013L	
27002KP	
27002K	27023VK
27002L	

Possible combinations: endoscope – cytology brush

Endoscope	Cytology brush
27000K	27023Y
27000L	
27001K	
27001L	
27002K	
27002L	
27003L	
27011K	
27011L	
27013L	
27002KP	

Possible combinations: light connection – fiber optic light cable

In combination with the product, light cables with a fiber diameter of 2.0 – 2.5 mm are recommended. The circle symbol on the product and the circle symbol on the light cable must match each other.

Circle symbols	Diameter of light cable
	4.8 – 5.0 mm
	3.0 – 3.5 mm
	2.0 – 2.5 mm

Combination: working channel – irrigation

Irrigation tubes can be connected via the LUER lock connection.














Possible combinations: eyepiece – camera head

Camera heads with a standardized KARL STORZ eyepiece connection can be connected to the eyepiece.

4.3 Ambient conditions

There are no special transport and storage conditions for this product.

4.4 Symbols on the packaging

Symbol	Meaning
	Manufacturer
	Date of manufacture
	Consult instructions for use
	Medical device
	Unique Device Identifier
	Article no.
	Batch code
	Serial number
	Number of products in the product packaging
	Unsterile
	Fragile, handle with care
	Federal (USA) law restricts this device to sale by or on the order of a physician.
	CE mark. The CE mark is a label required for certain marketable industrial products in accordance with EU law. It consists of the CE logo and, where applicable, the four-figure identification number of the responsible testing center (when the product belongs to a certain risk category).

5 Preparation

5.1 Unpacking the product

1. Carefully remove the product and accessories from the packaging.
2. Check the delivery for missing items and evidence of shipping damage.
3. In the case of damage, hidden defects, and short deliveries, document their nature and extent and contact the manufacturer or supplier immediately.

5.2 Assembling the product

1. Hold the instrument inlet on the two flat sides and push it onto the mount until the instrument inlet clicks into place.



2. If the instrument inlet does not click into place, push back the sleeve of the instrument inlet.
3. Insert the cock plug into the irrigation ports.



4. Place the spring cap on the thread of the cock plugs.
5. Tighten the spring cap slightly.
6. Fit the plug-in base (495F) onto the light inlet socket.
7. Screw the screw base (495G) onto the light inlet socket.

6 Disassembly

6.1 Disassembling the product

1. Unscrew the screw base from the light inlet socket.
2. Remove the plug-in base from the light inlet socket.
3. Unscrew the spring cap from the cock plug.
4. Remove the cock plug from the irrigation ports.



5. Push back the sleeve of the instrument inlet.
6. Remove the instrument inlet from the mount.



7 Maintenance, servicing, repairs, and disposal

7.1 Repairing the product

Repair work may only be performed by KARL STORZ or by a company authorized by KARL STORZ. The interventions described in this instruction manual are exempt from this rule.

- ▶ Please contact your local KARL STORZ subsidiary or authorized dealer (see the list of subsidiaries).

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, products must first be decontaminated. KARL STORZ reserves the right to send back contaminated products.

7.2 Disposing of the product

For disposal, the following measures are necessary:

1. Decontaminate the products prior to disposal.
2. Country-specific national laws and regulations must be observed.

8 Accessories and spare parts

8.1 Accessories

Item	Order no.
Fiber optic light cable, length 230 mm, diameter 2.5 cm	495NTA

8.2 Spare parts

Item	Order no.
Plug-in base, 9 mm in diameter, for Wolf fiberoptic light cables	495F
Screw base for KARL STORZ fiberoptic light cables and Olympus Corporation	495G
Seal for instrument inlets, pack of 10, single use recommended	27550N
Flow control stopcock	27504
Spring cap	6985691
Cock plug	8458190
Cleaning adaptor for instrument inlets	27001RA

9 Subsidiaries

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany
Postfach 230, 78503 Tuttlingen/Germany
Phone: +49 7461 708-0, Fax: +49 7461 708-105
Email: info@karlstorz.com

KARL STORZ Endoskope Berlin GmbH

Scharnhorststr. 3, 10115 Berlin/Germany
Phone: +49 30 3069090, Fax: +49 30 30 19452

KARL STORZ Endoscopy Canada Ltd.

7171 Millcreek Drive, Mississauga, Ontario L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free (Canada only) Phone: 1-800-268-4880, Fax: 1-800-482-4198
(Canada only)
Email: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.

2151 East Grand Avenue, El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free (USA only) Phone: 800 421-0837, Fax: 800 321-1304 (USA only)
Email: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.

1 South Los Carneros Road, Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
Email: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.

815 N. W. 57th Avenue, Suite 480, Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
Email: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.

Edificio Atlantic, Oficina 3G, Calle D e/ 1ra y 3ra, 10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
Email: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.

Av. Ejercito Nacional No. 453 Piso 2, Colonia Granada, Alcaldia Miguel Hidalgo, C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
Email: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.

Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi, CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
Email: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.

Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez, Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
Email: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS

Stamveien1, 1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
Email: post@karlstorz.no

KARL STORZ Endoskop Sverige AB

Storsåtragränd 14, 127 39 Skårholmen, Sweden
Phone: +46 8 505 648 00
Email: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY

Taivaltie 5, 01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
Email: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG

Representation Office
Kęstučio st. 59 / Lenktoji st. 27, 08124 Vilnius, Lithuania
Phone: +370 5 272 0448, Mobile: +370 685 67 000
Email: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S

Skovlytoften 33, 2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
Email: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.

415 Perth Avenue, Slough, Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
Email: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.

Displayweg 2, 3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
Email: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.

Phone: +31 (0)33 4545890
Email: info-be@karlstorz.com

KARL STORZ Endoscopie France S. A. S.

12, rue Georges Guynemer, Quartier de l'Europe, 78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
Email: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH

Landstraßer Hauptstr. 148/1/G1, 1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
Email: storz-austria@karlstorz.com

KARL STORZ Endoscopia Ibérica S. A.

Parque Empresarial San Fernando, Edificio Munich – Planta Baja, 28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
Email: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.

Via dell'Artigianato, 3, 37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
Email: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.

Capraška 6, 10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
Email: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.

Cesta v Gorice 34b, 1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
Email: pisarna@karlstorz.si

KARL STORZ Polska Sp. z o.o.

ul. Bojkowska 47, 44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
Email: info-pl@karlstorz.com

KARL STORZ Endoszkóp Magyarország Kft.

Toberek utca 2. fsz. 17/b, HU-1112 Budapest, Hungary
Phone: +36 195 096 31, Fax: +36 195 096 31
Email: info-hu@karlstorz.com

KARL STORZ Endoscopia Romania srl

Str. Prof. Dr. Anton Colorian, nr. 74, Sector 4, 041393 Bukarest, Romania
Phone: +40 (0)31 4250800, Fax: +40 (0)31 4250801
Email: info-ro@karlstorz.com

KARL STORZ Endoskope Greece M.E.P.E.*

Patriarhou Grigoriou E' 34, 54248 Thessaloniki, Greece
Phone: +30 2310 304868, Fax: +30 2310 304862
Email: info-gr@karlstorz.com

*Repair & Service Subsidiary

KARL STORZ Industrial**

Gedik Is Merkezi B Blok, Kat 5, D 38-39, Bagdat Cad. No: 162, Maltepe
Istanbul, Turkey

Phone: +90 216 442 9500, Fax: +90 216 442 9030

**Sales for Industrial Endoscopy

000 KARL STORZ Endoscopy – WOSTOK

Derbenyevskaya nab. 7, building 4, 115114 Moscow, Russia

Phone: +7 495 983 02 40, Fax: +7 495 983 02 41

Email: info-ru@karlstorz.com

TOV LLC KARL STORZ Ukraine

Avenue Geroyiv Stalingrada Str. 2D, office 717 Kyiv, 04210/Ukraine

Phone: +38 095 000-895-0, +38-097-000-895-0, +38 073 000-895-0

Email: marketing@karlstorz.com.ua

KARL STORZ SE & Co. KG Representation Office

Sabit Orudschow 1184, apt. 23, 1025 Baku, Azerbaijan

Phone: +99 450 613 30 60

Email: info-az@karlstorz.com

KARL STORZ ENDOSKOPE – East Mediterranean and Gulf (Offshore) S.A.L.

Spark Tower 1st floor Charles Helou St., Horch Tabet – Sin El Fil, Beirut,
Lebanon

Phone: +961 1 501105, Fax: +961 1 501950

Email: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.

P.O. 6061, Roggebaai, 8012 Cape Town, South Africa

Phone: +27 21 417 2600, Fax: +27 21 421 5103

Email: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan

Saryarka, 6, BC "Arman", off. 910, 010000 Astana, Republic of Kazakhstan

Phone: +7 7172 552-549, 552-788, Fax: -444

Email: info@karlstorz.kz

KARL STORZ ENDOSKOPE East Mediterranean & Gulf (branch)

Building West Side 7A – Unit 7WA – 3008, Dubai Airport Free Zone, P.O. Box
54983, Dubai - United Arab Emirates

Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282

Service Hotline: +971 (0)4 3415882

Email: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited

11th Floor, Dr. Gopal Das Bhawan, 28, Barakhamba Road, New Delhi
110001, India

Phone: +91 11 4374 3000, Fax: +91 11 4374 3010

Email: corporate@karlstorz.in

KARL STORZ SE & CO. KG

Interchange 21 Tower, Level 33, 399 Sukhumvit Road, North Klongtoey,
Wattana, 10110 Bangkok, Thailand

Phone: +84 28 3823 8000, Fax: +84 28 3823 8039

Email: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG

Resident Representative Office

14th Floor, MPlaza Saigon, 39 Le Duan, District 1, Ho Chi Minh City, Vietnam

Phone: +84 28 3823 8000, Fax: +84 28 3823 8039

Email: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.

Room 2503-05, 25F AXA Tower, Landmark East, No. 100 How Ming Street,
Kwun Tong, Kowloon, Hong Kong, People's Republic of China

Phone: +852 28 65 2411, Fax: +852 28 65 4114

Email: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch

Room 1805-1807, Building B, 18F Beijing IFC, No. 8, Jianguomenwai Street,
Chaoyang District, 100022, Beijing, People's Republic of China

Phone: +86 10 5638188, Fax: +86 10 5638199

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch

Room 701A Building 5 & Room 501 Building 7, No. 3000 Longdong Avenue,
Pilot Free Trade Zone, 201203, Shanghai, People's Republic of China

Phone: +86 21 60339888, Fax: +86 21 60339808

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch

Room 803-805, 8F Jin Jiang International Building, No. 1 West Linjiang
Road, Wuhou District, 6100414, Chengdu, People's Republic of China

Phone: +86 28 86587977, Fax: +86 28 86587975

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch

Room 2001-2005, 20F N-MEDIA International Center, No. 167 Youth Avenue,
Shenhe District, 110014, Shenyang, People's Republic of China

Phone: +86 24 23181118, Fax: +86 24 23181119

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch

Room 02B & 03 & 04A, 35F Teem Tower, No. 208 Tianhe Road, Tianhe
District, 510620, Guangzhou, People's Republic of China

Phone: +86 20 87321281, Fax: +86 20 87321286

Email: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.

No. 8 Commonwealth Lane #03-02, Singapore 149555, Singapore

Phone: +65 69229150, Fax: +65 69229155

Email: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd

No. 8 Commonwealth Lane #03-02, Singapore 149555, Singapore

Phone: +65 69229150, Fax: +65 69229155

Email: infoasia@karlstorz.com

KARL STORZ SE & Co. KG Representative Office Indonesia

Sinarmas MSIG Tower Level 37, Jl. Jend. Sudirman No. Kav. 21, Jakarta
Selatan

DKI Jakarta 12920

Email: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.

9F Hyowon-Building, 97, Jungdae-ro, Songpa-gu, 05719 Seoul, Korea

Phone: +82-70-4350-7474, Fax: +82-70-8277-3299

Email: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.

12F, No. 192, Sec. 2, Chung Hsin Rd., Sindian District, New Taipei City,
Taiwan

Phone: +886 933 014 160, Fax: +886 2 8672 6399

Email: info-tw@karlstorz.com

KARL STORZ SE & Co. KG Representative Office Philippines

1901 Picadilly Star Bldg., 4th Avenue, BGC, Taguig City 1636, Philippines

Phone: +63 2 317 45 00, Fax: +63 2 317 45 11

Email: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.

Stage Bldg. 8F, 2-7-2 Fujimi, Chiyoda-ku, Tokyo 102-0071, Japan

Phone: +81 3 6380-8622, Fax: +81 3 6380-8633

Email: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd .

68 Waterloo Road, Macquarie Park NSW 2113, P O Box 50 Lane Cove NSW
1595, Australia

Phone: +61 (0)2 9490 6700, Fax: +61 (0)2 9420 0695

Toll free: 1800 996 562 (Australia only)

Email: info@karlstorz.au

www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE



WIX212 • EN • V1.0 • 10-2020 • IFU • CE-MDR



KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Phone: +49 7461 708-0
Fax: +49 7461 708-105
E-mail: info@karlstorz.com
www.karlstorz.com

