

# Everything at a glance

SCENARA® .doorsign



## Always up-to-date with the digital OR door sign

The digital OR door sign from KARL STORZ offers detailed and automatically updated information about the room and procedure status as well as cleaning times, which is supported by an additional LED visualization of the OR status. It allows an overview of activities in the operating room without disrupting the OR workflow. Furthermore, the easy-to-use touchscreen provides a view of the entire daily schedule of the OR.




- 10" screen
- 255 x 181 x 25 mm
- VESA-compatible
- WiFi/LAN
- 12 V/PoE



## Everything at a glance

- Room number and room name
- Date and time
- Colored LED visualization of the OR status: free (light green), scheduled (light blue), occupied (blue), cleaning in progress (yellow)
- Visualization of risk information (red)
- Procedure details: Procedure type, date and scheduled time, patient status, progress indicator and additional information/comments
- Visualization of delays to improve time management and work processes
- Patient details (can be anonymized for data safety reasons): Name, date of birth, gender, infection risks
- Surgery team: Name and role
- Further room comments (e.g., maintenance)
- Current procedure status
- Display of cleaning times and start/stop triggers for cleaning management system



## SCENARA® – The modular application platform for OR processes

The digital door sign is part of the **SCENARA® portfolio**. As a modular and expandable software application platform, SCENARA® provides numerous applications for optimizing processes in and around the operating room.

The software complements existing solutions and expands and enhances their functionalities in order to simplify and optimize work processes. Interfaces to state-of-the-art hospital information and archiving systems allow seamless integration into the existing hospital infrastructure.

With SCENARA®, data security is ensured through an integrated user management, encryption in accordance with the latest security standards and a backup strategy.

More than  
**75**  
Years

*Shaping the Future  
of Endoscopy with you*